




Behind Every Smile is a Story

● QUALITY - ЯКІСТЬ ● PRECISION - ТОЧНІСТЬ ● ACCURACY - ПРАВИЛЬНІСТЬ ● DURABILITY - ЗНОСОСОТІЙКІСТЬ



- Super-sharp diamond-honed cutters
- Smooth ultra polished and contoured edges
- brilliant corrosion resistant stainless steel finish
- Superior T-12 stainless steel tool inserts
- sliding lock joint for balance & precision alignment

- Super-sharp diamond-honed cutters
- Smooth ultra polished and contoured edges
- brilliant corrosion resistant stainless steel finish
- Superior T-12 stainless steel tool inserts
- sliding lock joint for balance & precision alignment

STYLE
A

STYLE
B

ORTHODONTIC INSTRUMENTS

■ Инструменты ортодонтические
■ Instrumentos Ortodonticos
■ Orthodontische Instrumente

■ Инструменти ортодонтичні
■ Instruments Orthodontiques
■ ادوات تقويم الأسنان



786-1277
Aluminium Height
Gauge 0.22



786-1278
Aluminium Height
Gauge 0.18



786-1279
Bone Bracket Gauge



786-1280
Bone Thickness Gauge
40mm, 180mm



786-1281 OSTEOMETER
Measuring Forceps for Jaw
Ridge Width measurement, scale
on both sides, rustree, autoclavie



786-1282
Direct bonding
Bracket holders,
135mm



786-1283
Deluxe Bracket
holder Tweezers



786-1284
Direct bonding
Bracket holders,
136mm



786-1285
Direct Bond
Bracket holders,
Tweezers,
Small Tip



786-1286
Direct Bonding
Bracket holders,
140mm



786-1287
Direct Bonding
Bracket holders,
135mm



786-1288
Direct Bonding
Bracket holders,
140mm

■ Инструменты ортодонтические
■ Instrumentos Ortodonticos
■ Orthodontische Instrumente

■ Инструменти ортодонтични
■ Instruments Orthodontiques
■ ادوات تقويم الأسنان



786-1289
Band seater auto cleavable



786-1290
Ortho light scaler double ended



786-1291
Ortho light scaler double ended



786-1292
Ortho light scaler double ended



786-1293
Twiel-on elastic module instrument



786-1294
Twiel-on elastic module instrument



786-1295
Schure Band Seater



786-1296
Ortho light scaler



786-1297
Ortho light scaler



786-1298
Ortho light scaler heavy



786-1299
Ortho instrument



786-1300
Ortho instrument



786-1301
Ortho light Band Pusher Scaler


ORTHODONTIC INSTRUMENTS

■ Инструменты ортодонтические
■ Instrumentos Ortodonticos
■ Orthodontische Instrumente

■ Инструменти ортодонтични
■ Instruments Orthodontiques
■ ادوات تقويم الاسنان



786-1302
Ortho Light double ended
ligature director



786-1303
Ortho Light double ended
ligature scaler



786-1304
Ortho Light double ended ligature scaler director tivate coil spring
are remove modules with one instruments



786-1305
Ortho band pusher scaler




786-1306
Ortho Light Band Pusher scaler




786-1307
Ortho Ligature director single ended



786-1308
Ortho Light double ended ligature director scaler



786-1309
Ortho Light double ended ligature scaler director
tivate coil spring are remove modules with one instruments



786-1310
Ortho Light double ended ligature director



786-1311
Ortho Light ligature director scaler



786-1312
Ortho Light ligature director double ended



786-1313
Double ended cement spatula



786-1314
Ortho Light scaler single ended

■ Щипцы для ортодонтии и протезирования
■ Pincers pour l'orthodontie et la prothèse
■ Zanger für Kieferorthopädie und Prothetik


■ Щипці для ортодонції та протезування
■ Alicates para ortodoncia y prótesis
■ كماشة لتقويم الأسنان بروسشيتيس

BOX JOINT PLIERS

786-1315
Universal Dolphin Style
Wire Bending & cutting pliers



786-1321
Bracket Removing Plier curved
Bracket Removing curved for
easy removal of brackets in the
posterior region. Length 130mm



786-1316 Band Removing Plier
The curved blade of the band
removing plier can strip easily
and gently remove directly to
the tooth, Length 130mm



786-1322
Young Wire Bending &
Loop Forming Pliers



786-1317
Tweed Loop
forming Pliers



786-1323
Wire Bending &
Loop Forming Pliers



786-1318 How Plier Curved
Universal usable tool for
the multiband technology for
wires up to .028"
Length 130mm



786-1324
Adams, small wire bending
and loop forming Pliers



786-1319 How Plier
Universal usable tool for
the multiband technology for
wires up to .028"
Length 135mm




786-1325
Bracket Removing Plier Round Box
Bracket removing plier for comfortable
and quick remover of the brackets.
Length 135mm



786-1320
Nance Loop
forming Pliers



786-1326
Weingart Plier medium
Weingart plier that are ideal due to
their curved tip for setting sheets and
is capable of bending the ends
Length 145mm



PLIERS FOR ORTHODONTICS & PROSTHETICS




- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zanger für Kieferorthopadie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كاشاة لتقويم الأسنان بروسشيتيس

BOX JOINT PLIERS


786-1327

Adams Plier Round Box
Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak Length 125mm



786-1333 Adams Plier

Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak Length 125mm



786-1328

Falcon Beak Contouring Pliers




786-1334 Tweed

Loop Forming Pliers




786-1329

Wires up to 0.020" / 0.7mm




786-1335 Bird Beak Plier Standard

Bird Beak Wire bending plier without notches, which is suitable as a universal utility too Length 125mm



786-1330

Johnson Contouring Pliers



786-1336

Mini Ligature cutters



786-1331 Aderer

Large wire bending and loop forming pliers




786-1337 Small Fluch-cut distal and cutting Pliers

with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



786-1332 Tweed

Angle rectangular arch forming & Torquing pliers



786-1338 Large Fluch-cut distal and cutting Pliers

with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



■ Щипцы для ортодонтии и протезирования
■ Pincers pour l'orthodontie et la prothèse
■ Zanger für Kieferorthopädie und Prothetik

■ Щипці для ортодонції та протезування
■ Alicates para ortodoncia y prótesis
■ كماشة لتقويم الأسنان بروتستيتيس

BOX JOINT PLIERS

786-1339
Tweed / Angle
1.78mm (.070") wide tips arch
forming and torquing pliers



786-1345 Large
Contouring Plier



786-1340
Tweed / Angle
1.78mm (.070") wide tips arch
forming and torquing pliers



786-1346
Bands & Crowns Plier
In order to skirt



786-1341
Lingual / Palatal arch
forming pliers



786-1347
Ligature Forming Pliers



786-1342
Aderer
Facebow loop forming plier



786-1348
For double back and triple back
bands in 0.30" / 0.7mm & 0.36
/0.91mm lingual arch wires



786-1343
Waldsavhs wire bending
and cutting pliers





786-1349
**Bird Beak Wire Bending
& cutting pliers**



786-1344 Aderer Plier
For the manufacturing of appliances and
activation of wires for wires up to .030" fine
tip and up to .036" mean tip
Length 130mm



786-1350
Weingart Plier Medium
Weingart plier that are ideal due to
their curved tip for setting sheets and
is capable of bending the ends.
Length 145mm



PLIERS FOR ORTHODONTICS & PROSTHETICS



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zanger für Kieferorthopadie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- کماشة لتقويم الأسنان بروسشیتیس

L-KEY JOINT PLIERS



786-1351
Distal end cutter with safety hold mini head, long handles. fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.40 mm to 0.64 mm (0.18" x .025")
- ⊗ 0.38 mm to 0.53 mm (0.155" to .0215")
- ⊗ Max. 0.53 mm to 0.64 mm (0.21" x .025")



786-1352
Flush-cut distal end cutter with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.46 mm to 0.64 mm (0.18" x .025")



786-1353
Soft wire and ligature cutter Mini fitted with un-breakable corrosion resistant and durable stainless steel tips

- ☞ .25 mm to .38 mm (.010" to .015")



786-1354
Soft wire and ligature cutter, standard fitted with un-breakable corrosion resistant and durable stainless steel tips

- ☞ .25 mm to .38 mm (.010" to .015")



786-1355
Soft wire and ligature cutter 15° angled fitted with un-breakable corrosion resistant and durable stainless steel tips

- ☞ .25 mm to .38 mm (.010" to .015")



786-1356
Soft wire and ligature cutter 45° angled fitted with un-breakable corrosion resistant and durable stainless steel tips

- ☞ .25 mm to .38 mm (.010" to .015")



786-1357
Hard wire cutter, straight fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1358
Hard wire cutter, 15° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1359
Hard wire cutter, 15° angled Non-autocleavable

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.56 mm to 0.71 mm (0.22" x .028")

■ Щипцы для ортодонтии и протезирования
■ Pincers pour l'orthodontie et la prothèse
■ Zanger für Kieferorthopädie und Prothetik

■ Щипці для ортодонції та протезування
■ Alicates para ortodoncia y prótesis
■ کماشة لتقويم الأسنان بروسشیتیس

L-KEY JOINT PLIERS



786-1360
Band removing pliers



786-1365
Weingart-mini
Utility Pliers




● 0.30 mm to 0.50 mm (0.12" to .020")
■ Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1361
Bracket removing
Pliers angled



786-1366
Weingart-mini
Utility Pliers



● 0.30 mm to 0.50 mm (0.12" to .020")
■ Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1362
Bracket removing pliers straight



786-1367
How Straight
Utility Pliers




● 0.30 mm to 0.50 mm (0.12" to .020")
■ Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1363
Bracket removing pliers with plastic tip



786-1368
Angled
Utility Pliers



● 0.30 mm to 0.50 mm (0.12" to .020")
■ Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1364
Adhesive removing pliers



786-1369
Wire bending and loop forming pliers
Angled (Bird Beak)



● 0.30 mm to 0.50 mm (0.12" to .020")
■ Max. 0.56 mm to 0.71 mm (0.22" x .028")



PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zanger für Kieferorthopadie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- کماشة لتقويم الأسنان بروسشیتیس

L-KEY JOINT PLIERS



786-1370
Multi-purpose pliers, bird beak
De la rosa and cutter in one

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.46 mm to 0.56 mm (0.18" x .022")



786-1375
Stop pliers V-Bebd - 1mm

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.43 mm 0.56 mm (0.17" x .028")



786-1371
Cinch Back / Tie back
Distal end bending pliers

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1376
Lingual / palatal arch forming pliers

- 0.76 mm 0.30" / 0.91 mm (.036")




786-1372
Hammer head
Distal end bending pliers

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1377
Torquing pliers set



786-1373
Nance Loop forming pliers

- 0.36 mm to 0.71 mm (0.14" to .028")
- Max. 0.56 mm to 0.70 mm (0.22" x .028")




786-1378
Tweed / angled
Rectangular arch forming and
torquing pliers

- 0.30 mm to 0.71 mm (0.12" to .028")
- Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1374
Detailing step pliers / Bayonet
bending pliers - 1mm

- 0.30 mm to 0.50 mm (0.12" to .020")
- Max. 0.56 mm to 0.71 mm (0.22" x .028")



786-1379
Surgical ball hook crimping pliers
angled




786-1380
Surgical ball hook crimping pliers
straight



■ Щипцы для ортодонтии и протезирования
■ Pincès pour l'orthodontie et la prothèse
■ Zanger für Kieferorthopädie und Prothetik

■ Щипці для ортодонції та протезування
■ Alicates para ortodoncia y prótesis
■ كماشة لتقويم الأسنان بروتستيتيس


L-KEY JOINT PLIERS





786-1381
Aderer Clamp three jwa wire bending plier for contouring and bending wire up to 30" / 0.76mm




786-1387
Falcon beak contouring pliers





786-1382
Bird beak plier for bending round square and rectangle wire up to 30" / .76mm




786-1388
For placement and removal of arch wires, pins and other auxiliaries





786-1383
Bird beak plier for light wire bending round square and rectangular up to .016" / 0.41mm



786-1389
For placement and removal of arch wires, pins and other auxiliaries in angle design





786-1384
Double Rounded jaw bending pliers hooks and wire up to .020" / 0.51mm




786-1390
Contouring Plier




786-1385
Adams pliers to hooks and bending wire up to .030" / 0.76mm




786-1391
Bending plier, in order to skirt molaras bands, with grooved Recommended wire 016". 018" and 022" .041mm, 0.46mm and 0.56mm



786-1392
Long forming in four length 3,4,5, 7 6mm wire capacity 0.22 . 56mm



786-1386
Johnson pliers for crown band contouring




786-1393
Arch bending plier design for convenient handling of square or rectangular wire, 020" / 0.51mm






- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zanger für Kieferorthopadie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- کماشة لتقويم الأسنان بروتھتیس



LAP & L-KEY JOINT PLIERS



786-1394
In order to skirt bands and crowns.




786-1395
Wires up to .020" / 0.7mm




786-1396
Round and concave plier for loop forming and bending wire up to .20" / .051mm




786-1397
Bending for round square or rectangular wire up to 0.30" / 0.76mm



786-1398
for double back and triple back bands in 0.30" / .76mm & .036" / 91mm lingual arch wires



786-1399
Detailing step plier double sided offset beak permits bayonet bends. Creates right hand or left hand step-up or stepdown at 1/2", 3/4" or 1mm increments.



786-1400
"V" bend stop special design makes 1mm "V" bends to shorten archwires or provide a positive stop Excellent for placing stop in nitanium wires For round of rectangular wires up to .025"



786-1401
Angle how plier efficient utility plier with serrated beak for a firm grip 45 angled tips to access hand or reach areas.



786-1402
Weingardt plier serrated tips hold wires firmly at a convenient working angle tapered beak fit easily between brackets and hard to reach areas, Our most popular utility plier



786-1415
Pin 8 Ligature Cutter TC




786-1416
Mini Pin 8 Ligature Cutter TC



786-1417
Standard Ligature Cutter TC



786-1418
V Band Stop Plier TC.



786-1419
Angled Direct Bond Bracket Remover TC



786-1420
Distal end Cutter Safety hold Plier
(Ideal for Cutting Wire) TC



786-1421
Hook Crimping Plier Angled
Straight TC



786-1422
3 Prong Clasp Adjustable
Plier TC



786-1423
Utility Plier Weingart TC



786-1424
Diret Bracket Remover
TC



786-1425
Straight Direct Bond Bracket
8 Remover TC



786-1426
Pin Holding Plier TC



786-1427
Hard wire cutter



786-1428
Elastic Separating Plier
Round Handle



786-1429
Elastic Separating Plier



786-1430
786-1431
Smaha
Elastomeric
ligature
placing forceps
140mm



786-1432
Elastomeric
ligature placing
forceps with
groove
Falcon 140mm



786-1433
Elastomeric
ligature placing
forceps with
groove
Falcon 140mm

